

PLHX1002

Respiratory endoscope training model

01

Respiratory endoscope
intervention training

Main components

Mobile head, nasal cavity, larynx (epiglottis and glottis), main airway, bronchi (second to fifth-level branches)

Product function

It is applicable to the performance test, functional demonstration and training of interventional devices, such as laryngoscope, bronchoscope/ soft endoscope and airway biopsy devices

Items supported

- Bronchoscopy operation training
- Anatomical structure identification under endoscope
- Bronchoscopic device operation training

Main material

SLD silica gel

