

PLMN1001

Urological endoscope training model

03

Uroscope intervention training

Main components


Male reproductive organs (penis, testis), urethra, prostate, urinary bladder, ureter, kidneys (calyces, renal pelvis), simulated calculi, etc.

Product function

It is applicable to the performance test, functional demonstration and training of intervention devices, such as cystoscope, ureteroscope, and guidewire/ catheter

Items supported

- Pyelolithotomy and ureterolithotomy training;
- Ureteroplasty (dilatation and stenting) training;
- Urethroplasty (dilatation and stenting) training;
- Male urinary bladder irrigation training
- Male catheterization training;
- Transurethral resection of prostate training;
- Calculi search and extraction training under nephroscope;
- Simulated calculi (optional components) can be placed in the renal pelvis according to the training demands to practice the flexible ureteroscopic holmium laser lithotripsy



Main material

SLD silica gel



Accessories: Renal calculus



Note: kidneys can be customized according to demands